

Preparation of chemically competent *Escherichia coli* cells

Materials

0.5 or 1.5-ml microfuge tubes
50-ml Falcon tubes

Chemicals

DMSO

Procedure

1. Inoculate 5 ml LB medium with the appropriate antibiotic(s) with the *E. coli* strain of which you want to make competent cells and incubate overnight at 37°C.
When preparing DH5 α competent cells it is better to use SOB medium instead of LB.
2. Use the overnight culture to inoculate 500 ml LB medium and incubate at 30°C until the absorbance at 600 nm is between 0.4-0.6.
Optional: add 2.5 ml 2M MgCl₂ to the medium (to a final concentration of 10 mM) at the start of the cultivation.
In the original paper (Inoue et al, 1990) the culture was incubated at 18°C but similar results were obtained when the cultures were incubated at higher temperatures (25-30°C). SOB medium for DH5 α competent cells.
3. Chill the culture for at least 10 min on ice.
In the following steps the cell suspension should be kept on ice as much as possible.
4. Spin the cell suspension for 10 min at 6000 rpm (Sorvall GSA rotor) or 4000 rpm when harvested in 50-ml Falcon tubes.
5. Gently resuspend the pellet in 100 ml ice-cold CC buffer in 50-ml Falcon tubes.
Resuspend with a 10-ml serological pipette and avoid introducing bubbles.
6. Incubate the cell suspension on ice for at least 10 min.
7. Spin for 10 min at 4000 rpm at 4°C.
8. Gently resuspend the pellet in 18.6 ml ice-cold CC buffer and add 1.4 ml DMSO.
9. Incubate the cell suspension on ice for at least 10 min.
10. Distribute the cell suspension in 100-200 μ l aliquots in 0.5 or 1.5-ml microfuge tubes.
11. Flash freeze the cell suspension in liquid nitrogen and store the tubes at -80°C.
At -80°C the cell will be competent for at least 6 months.

CC buffer

| | | |
|--------|--------------------------------------|-----------|
| 10 mM | Hepes | 2.38 g/L |
| 15 mM | CaCl ₂ | 2.21 g/L |
| 55 mM | MnCl ₂ ·4H ₂ O | 10.89 g/L |
| 250 mM | KCl | 18.64 g/L |

Dissolve all components except MnCl₂ and adjust the pH to 6.7 with KOH. Then add the MnCl₂ and filter sterilize the solution over a 0.22 μ m filter.

In some protocols Pipes is used instead of Hepes.

2 M MgCl₂

| | |
|---------------|--------------------------------------|
| 40.66 g/100mL | MgCl ₂ .6H ₂ O |
|---------------|--------------------------------------|

Dissolve 40.66 g MgCl₂.6H₂O in 74 ml water. Sterilize by autoclaving and store at room temperature.

LB medium

| | |
|--------|---------------|
| 10 g/L | tryptone |
| 5 g/L | yeast extract |
| 10 g/L | NaCl |

Sterilize by autoclaving and store at room temperature.

SOB medium

| | |
|-----------|---------------|
| 20 g/L | tryptone |
| 5 g/L | yeast extract |
| 0.5 g/L | NaCl |
| 0.186 g/L | KCl |

Adjust the pH to 7.0 with NaOH.
Sterilize by autoclaving and store at room temperature.

Reference

Inoue *et al.* (1990) *Gene* 96, 23-28.